

## Your SureMed or Qube Mini Starter Kit



Your SureMed or Qube Mini kit includes: 26 blisters & cards, filling foam & sealing roller.

## Tips


- Using the reverse side day and time of day as a guide fill the appropriate blister with the correct medication.
- This blister packaging is not child resistant, please store beyond the reach of children.
- Keep unused cards in a clean, dry and air condition location.
- If there is a slight tug while peeling off the pressure sensitive liner try peeling from a different corner.

### 1. Remove Liner



Remove the pressure sensitive liner revealing the sealing surface.  
**NOTE: Qube Card has a pressure liner on the middle panel only.**

### 2. Filling Foam



Place the card on the filling foam template (dry glue side up).

### 3. Insert Blisters



Insert the blister sheet into the card starting with the top row. Lower the rest of the sheet into the remaining holes.

### 4. Filling Pill #1



Say you take 7 medications a day and one pill twice a day. Pill #1 you take in the morning. Place one pill in each of the morning blisters, Sunday thru Saturday.

### Filling Pill #1



You also take the same (pill #1) each evening as well. Place one pill in each of the evening blisters, Sunday thru Saturday.

### Filling Pill #2



Pill #2 you take each morning. Place one pill in each of the morning blisters, Sunday thru Saturday.